



# Enhancing Implant Fixation by using Porous Plasma-Sprayed Titanium Coatings on Metal, Ceramics & PEEK



## TPS METAL

Plasma-sprayed titanium coatings have been applied to medical implant surfaces for several years. The coating process has proven to be a very successful means of producing highly adherent, micro-porous and biocompatible implant surfaces [1, 2, 3].


Porous titanium coatings have been plasma-sprayed onto orthopedic implant surfaces – for clinical applications – since the early 1980 [4]. As a result of their structure, the surfaces ensure that the implant is firmly anchored in the bone.



## TPS CERAMIC

An increasing number of implants and implant components are being manufactured from ceramic substrates; which are known to offer excellent biomechanical properties and high mechanical strength. Ceramic implants are properly classified as bioinert [5] and therefore do not foster cementless fixation. The ceramic implant surface must therefore be modified to promote osseointegration.

Porous TPS coatings are applied to the implant by means of a vacuum plasma spray process with which the surface of the ceramic implant is modified. The topography of the implant surface is thus optimized for bone apposition, providing bone cells with the ideal topography for ingrowth without reducing the fatigue strength of the implant body.




## TPS PEEK

In addition to coating metals and ceramics, we offer our customers a highly porous TPS coating for implant components based on temperature resistant polyethylene (PEEK). The TPS coating is applied onto the implants by using vacuum plasma spraying

technology, actively yielding a high degree of porosity and coating adhesion. Biocompatible surfaces, which are attractive to the bone, can also be applied to substrates fabricated from bioinert plastics. This in turn, causes osseointegration to occur.




Test criteria TPS METAL	Result
Color	Gray
Coating thickness (ASTM F 1854)	30-800 $\mu\text{m}$
Tensile strength (ASTM F 1147)	$\geq 22$ MPa
Porosity (ASTM F 1854)	20-40 %
Roughness $R_a$ (DIN EN ISO 3274) (DIN EN ISO 4287) (DIN EN ISO 4288)	3.5-80 $\mu\text{m}$
Biocompatibility (DIN EN ISO 10993-1)	Biocompatible
Titanium powder	In accordance with ASTM F 1580




Vertical micrograph of the TPS METAL coating

Test criteria TPS CERAMIC	Result
Color	Gray
Coating thickness (ASTM F 1854)	50-300 $\mu\text{m}$
Tensile strength (ASTM F 1147)	$\geq 22$ MPa
Porosity (ASTM F 1854)	20-40 %
Roughness $R_a$ (DIN EN ISO 3274) (DIN EN ISO 4287) (DIN EN ISO 4288)	30-60 $\mu\text{m}$
Biocompatibility (DIN EN ISO 10993-1)	Biocompatible
Titanium powder	In accordance with ASTM F 1580



SEM picture of the TPS CERAMIC coating

Test criteria TPS PEEK	Result
Color	Gray
Coating thickness (ASTM F 1854)	50-250 $\mu\text{m}$
Tensile strength (ASTM F 1147)	$\geq 22$ MPa
Porosity (ASTM F 1854)	20-60 %
Roughness $R_a$ (DIN EN ISO 3274) (DIN EN ISO 4287) (DIN EN ISO 4288)	20-40 $\mu\text{m}$
Biocompatibility (DIN EN ISO 10993-1)	Biocompatible
Titanium powder	In accordance with ASTM F 1580



Micro-section of the TPS PEEK coating

### References

1. Gruner H. (2003) Plasmaspritzen für die Medizintechnik. medrevat.ch/deutsch/plasmaspritzen.htm
  2. van Osten U.<sup>1</sup>, Salito A.<sup>2</sup>, Breme F.<sup>1</sup>, Aits M.<sup>1</sup>, Hufnagel K.<sup>1</sup> (1995) Maßschneidern poröser Titanschichten für die zementfreie Implantationstechnik. <sup>1</sup>GfE Metalle und Materialien GmbH, Nuremberg, Deutschland. <sup>2</sup>Sulzer Metco AG, Wohlen, Schweiz
  3. Lembke U. (2010) Implantatbeschichtung mit Vakuum-Plasma-Verfahren. Referat aus dem Fachforum „Implantate-Materialien, Beschichtungen und Anwendungen“. Ostbayrisches Technologie-Transfer-Institut e.V. Bayreuth
  4. Luedemann R. E. (1996) Mechanical Characterization of the Interaction between a plasma sprayed Ti alloy coating and Ti alloy and Co-Cr Substrates. 42nd Annual Meeting, Orthopaedic Research Society, February 19-22.
  5. Giulio Maccauro, Pierfrancesco Rossi Iommetti, Luca Raffaelli and Paolo Francesco Manicone (2011) Alumina and Zirconia Ceramic for Orthopaedic and Dental Devices, Biomaterials Applications for Nanomedicine, Prof. Rosario Pignatello (Ed.), ISBN: 978-953-307-661-4, InTech
- Devine D.M., Hahn J., Richards R.G., Gruner H., Wieling R., Pearce S.G. (2013) Coating of carbon fiber-reinforced polyetheretherketone implants with titanium to improve bone apposition. J Biomed Mater Res B Appl Biomater. 101(4): 591-598.
  - Walsh W.R. et al. (2015) Plasma-sprayed titanium coating to polyetheretherketone improves the bone-implant interface. The Spine Journal
  - Acklin Y.P. et al. (2016) Prospective clinical and radiostereometric analysis of the Fitmore short-stem total hip arthroplasty. Arch Orthop Trauma Surg. 136(2): 277-284.
  - Bouras et al. (2015) Good long term results following cementless TKA with a titanium plasma coating. Knee Surg Sports Traumatol Arthrosc. DOI 10.1007/s00167-015-3769-3.

We look forward to talking with you!



DOT GmbH  
A Member of Eifeler Holdings  
Charles-Darwin-Ring 1a  
18059 Rostock

Phone: +49 381 40 33 50  
Fax: +49 381 40 33 599  
info@dot-coating.de  
www.dot-coating.de

### DOT – coating specialist for orthopedic and dental implants

DOT is one of Europe's leading providers of medical coating solutions for orthopedic and dental implants and instruments and also their cleanroom packaging.

Our comprehensive supply chain concept makes us an ideal medical technology partner. Our activities help restore the health of patients worldwide and thus make a major contribution to the improvement of their quality of life.